















Cathejell[®] with Lidocaine

Product	Indication	Pack Size
<p>CATHEJELL[®] with Lidocaine</p> <ul style="list-style-type: none"> - Catheter lubricant gel with anesthetic <p><i>Contains: Lidocaine, 2g Chlorhexidinedihydrochloride 0.05g (accordion-type)</i></p>	<p>Instillation into urethra prior to introduction of catheters or other instruments.</p> 	<p>Individually Blister Pack, 25 tubes of 12.5gm each</p>

Sterile  gel in sterile  packaging. Steam  sterilized, therefore no residues of ethylene oxide in the product. Sterility and the disinfecting agent  Chlorhexidine help to avoid iatrogenic infections (important for indwelling catheters). The  collapsible tube with smooth anatomic nozzle and internal  break-off tip (patented, no danger of trauma) allows  soft instillation as well as aspiration of urine if necessary.  Watersoluble,  crystalclear.  Tissue adapted conductivity: ideal in combination with metal-shaft  instruments.